

NeuroAnalyzer report: Extract LDF Data

NeuroAnalyzer v0.3.0 · MATLAB 2026a · 2026-09-26T13:37:08+00:00 · window image: exportapp

Extract LDF Data

Extract LDF export (demo) in time window, crop (in demo)

? Help

2 Choose time range

Type Start (s), end (s) "Pick on plot" and then click the start and end on either plot

start (s)

20.000

end (s)

280.000

Pick on plot

Full range

3 Crop

Crop to range

Show

Cropped segment

Cropped 20.000-280.000 s
4 min 20.0 s, 260001 samples

4 Save

Save cropped data...



Croped [stim, LDF, Fz] as cropped LDF Processing

Save session...

Open session...

Report (PDF)...

Methods text...



✓ Session saved to Extract_LDFApp_session.mat (1 input file(s) with MTF, settings and results).

© Copyrights by Alejandro Suarez, Ph.D. - v0.3.0

Window: Extract LDF Data (ExtractLDFApp)
NeuroAnalyzer v0.3.0
MATLAB 26.1.0.3346908 (R2026a) Update 5
OS: Linux 6.17.0-1022-azure (GLNXA64)
Session created: 2026-09-26T13:37:08+00:00
Report created: 2026-09-26T13:37:08+00:00

NOTES
ExtractLDFApp walkthrough: demo data

INPUTS (1)
1. LDF export: demo_ldf_export.mat
path: /tmp/NeuroAnalyzerDemo/demo_ldf_export.mat
18,310,681 bytes, modified 2026-09-26T13:29:39
MD5 3b9e6147ced8d96db6949e2b1b947853

SETTINGS
range = [20 280]
view = Cropped segment
stimulusChannel = 6
ldfChannel = 8

KEY RESULTS
Recording: 300000 samples at 1000 Hz (5 min 0.0 s)
Crop: 20.000-280.000 s, 260001 samples
Cropped LDF: mean 127.7, SD 10.44