

NeuroAnalyzer report: Extract Ephys Data

NeuroAnalyzer v0.3.0 · MATLAB 2026a · 2026-09-25T17:59:28+00:00 · window image: exportapp

Extract Ephys Data

Local recording (TDT, Biotin, Open Ephys, NWB), extract LFP and MUA

Help

RAW: 24414.06 Hz · stim: 1017.25 Hz · 33.0 s
8 RAW ch · 2 stimuli

2 Choose channels

Stimulus channel

Ch 1

RAW channels

All

None

Ch 1
Ch 2
Ch 3
Ch 4
Ch 5
Ch 6
Ch 7

3 Process

Process LFP...

Process MUA...

Plot RAW

4 Save

Save LFP...

Save MUA...

Export NWB...

Save session...


Open session...

Report (PDF)...

Methods text...

LFP: 8 ch @ 1017.25 Hz

Processed LFP (1017.25 Hz)



✓ Session saved in ExtractEphysApp as session named '1 (input Biotin)' with NWB settings and results.

© Copyright by Alejandro Suarez, Ph.D. · v0.3.0

Window: Extract Ephys Data (ExtractEphysApp)

NeuroAnalyzer v0.3.0
MATLAB 26.1.0.3346908 (R2026a) Update 5
OS: Linux 6.17.0-1022-azure (GLNXA64)
Session created: 2026-09-25T17:59:28+00:00
Report created: 2026-09-25T17:59:28+00:00

NOTES

ExtractEphysApp walkthrough: demo data

INPUTS (1)

1. TDT tank recording: demo_tank
path: /tmp/NeuroAnalyzerDemo/demo_tank
21,458,904 bytes, modified 2026-09-25T17:55:03
MD5: b7a69c4c9f27f2332fba7c86b6b358ae

SETTINGS

format = tdt
stimChannel = 1
rawChannels = [1 2 3 4 5 6 7 8]
lfp.params.lowCutoff = 300
lfp.params.notch60 = false
lfp.params.downsample = true
lfp.params.downsampleRate = 1000
lfp.channels = [1 2 3 4 5 6 7 8]
lfp.stimChannel = 1
mua = []

KEY RESULTS

Recording demo_tank: 8 RAW channels at 24414.06 Hz
LFP: channels [1 2 3 4 5 6 7 8] at 1017.25 Hz, 30518 samples, RMS
[1.55e-05 1.59e-05 1.72e-05 1.82e-05 1.73e-05 1.59e-05
1.55e-05 1.54e-05] V